

Nextera® Library Validation and Cluster Density Optimization

Guidelines for generating high-quality data with Nextera library preparation kits.

Introduction

Each kit in the Nextera family of library preparation (DNA, XT DNA, Custom Enrichment, and Exome Enrichment) requires different considerations for final library validation and optimization of cluster density. Because each Nextera protocol is designed for a distinct application, it is important to tailor the library validation and clustering processes to each specific Nextera kit.

The Nextera “tagmentation” reaction simultaneously fragments and tags DNA with adapters. As with other enzymatic fragmentation protocols, the amount of input DNA can affect the outcome of the assay. Using too little DNA may result in libraries with shorter insert sizes, and conversely using too much DNA will result in libraries with larger insert sizes. This technical note describes best practices for optimizing the final libraries generated with each Nextera kit.

Starting Material

Nextera kits are ideal for high-throughput studies, providing fast library preparation methods with low input requirements. Nextera DNA and Enrichment protocols are optimized for 50 ng of high-quality genomic DNA, and the Nextera XT DNA kit requires 1 ng of small genome or amplicon DNA. To ensure that nucleic acids aside from double-stranded DNA do not interfere with quantitation, Qubit¹ or PicoGreen² assays are recommended to quantify the input genomic DNA. To ensure accuracy, Illumina recommends quantifying a concentrated solution of each sample and subsequently diluting to the input concentration specified in each library preparation user guide.

Organic contaminants (including ethanol) have been shown to interfere with the Nextera tagmentation reaction, and should be eliminated from the input DNA prior to quantitation using standard methods. EDTA in DNA elution or dilution buffers can also interfere with the tagmentation reaction, and should be avoided.


Nextera DNA Library Preparation

Library Validation

A wide range of final library sizes is expected with Nextera DNA protocol, and the yield and size profile will be highly dependent on the amount of input DNA. Typical libraries demonstrate a broad size distribution, with average fragment sizes ranging from 250 bp to 1.5 Kb in length. Figure 1 shows two examples of final library traces.

Illumina recommends using a BioAnalyzer³ to perform quality control of the library and determine library size. A Qubit or PicoGreen assay should be used to quantify the double-stranded final library. One of

Figure 1: Distribution of Final Library Sizes



The final library traces for Nextera DNA libraries were generated using a BioAnalyzer. These example size distributions are also applicable to Nextera XT libraries that have not been normalized. To optimize clustering concentrations based on average library size of 500 bp, the sample shown in Panel A requires the conversion factor 1 ng/μL = 3 nM. Clustering concentrations for the sample shown in Panel B should be converted using the equation 1 ng/μL = 1.5 nM, based on an average library size of 1 Kb.

three conversion factors for library concentration is recommended, depending on the average library size. Based on an average size of 500 bp, library concentrations generated by the quantification step can be converted to nM with the equation **1 ng/μL = 3 nM**. For libraries with average size ≥ 1 Kb, the conversion factor **1 ng/μL = 1.5 nM** is recommended. For shorter libraries (~250 bp), a conversion factor of **1 ng/μL = 6 nM** should be used (Table 1). Libraries with BioAnalyzer profiles similar to the examples shown in Figure 1 are expected to yield high-quality sequencing data, once appropriate clustering conditions have been established (see below).

